

EUGENAT

Liquid for root canal desobturation

Composition

Thymol
Isoamyl acetate
Tetrachlorethylene

Indications and properties

In cases of repeated endodontic treatment it is necessary to clean the root canal from the old filling material. To do this, use a special liquid to facilitate such manipulations. "Eugenat" is used for softening the pastes based on zinc oxide and eugenol.

Recommended use

Method for unsealing the root canals in one visit.

Proceed with the mechanical treatment of pulp chamber and the entrance to the root canals as widely as possible. Using a probe introduce one drop of "Eugenat" in a tooth cavity, and then into the root canals and leave it for 1 minute in order to soften a paste. Then begin processing with the endodontic instrument every time wetting its tip with "Eugenat".

After passing 1/2 of the root canal re-fill the canal with "Eugenat" liquid and leave for 1 minute with the aim of softening a paste in the deep part of the root canal. Then start treatment the canal with endodontic instruments every time wetting its tip with "Eugenat" as long as you reach the apex. After unsealing the root canal should be washed, degreased, and dried in order to prepare for subsequent filling.

Method for unsealing the root canals in two visits.

Proceed with the mechanical treatment of pulp chamber and the entrance to the root canals as widely as possible. Wet a cotton tampon with "Eugenat" and place it to the root canal, slightly compacted. Close the root canal with temporary filling until the next visit. After 1-3 days remove the temporary filling and proceed with mechanical unsealing the root canal from softened paste that in this time would be easy. After unsealing the root canal should be washed, degreased, and dried in order to prepare for subsequent filling.

Package contents and storage

Liquid is available in a glass bottle of 13 ml.
Keep in cool dry place away from direct sunlight.
Shelf life is 3 years.